

Honeywell Specialty Chemicals Seelze GmbH

Postfach 10 02 62, 30918 Seelze

Tel.: +49 (5137) 999-0, Fax: +49 (5137) 999-123

www.honeywell.com

Honeywell Specialty Chemicals Seelze GmbH · Postfach 10 02 62 · 30918 Seelze

Ihre Zeichen

Ihre Nachricht

Unsere Zeichen

Telefon: 05137 999-

Seelze, 26.04.2021

Hydranal CoE

353

HYDRANAL™ NEXTGEN FA Reagents

To Whomsoever It May Concern

This letter confirms that there is no need to develop new analysis methods when upgrading to HYDRANAL NEXTGEN FA reagents.

The suffix "FA", added to the name of the reagent, indicates that the reagent is "free of alcohols". Alcoholic solvents that are typically used in standard KF reagents have been replaced by acetonitrile in combination with confidential stabilizers and accelerators which guarantee an uncompromising performance of Hydranal reagent. Customers who wish to upgrade their Hydranal reagents to the corresponding Hydranal FA reagents can do so without changes in the titration method.

For more detailed information please do not hesitate to contact hydranal@honeywell.com.

Kind regards,

Honeywell Specialty Chemicals Seelze GmbH

Thomas Wendt
Supervisor HYDRANAL Center of Excellence



Sabine Mayer
Quality Manager